

SecureConnect® needle-free access

Needle-free access optimised for safety

Septum surface easily disinfected

Smooth, gapless surface

Septum seal reduces risk of biofilm formation

Housing tight between septum

Fluid pathway facilitates optimum flushing

Straight and direct

Zero dead space

Designed to prevent harbouring of contaminants and biofilm

Straight through flow


Maximises flow rates

No moving parts reduces risk of HA-BSI* contamination

No internal mechanism

Clear visibility promotes complete and adequate flush

Transparent polycarbonate and BPA**-free copolyester



* Hospital Acquired Blood Stream Infection

**Bisphenol A (BPA)