

TIVA TECHNICAL SPECIFICATIONS



Colour-coded Anti-Syphon valves and clamps

- Prevent accidental free-flow
- Prevent back-tracking of fluids and drugs
- Clear identification of lines

PE lined tubing

- All drug lines are PE-lined tubing
- Prevents drug binding
- No leaching of phthalates into fluid pathway

Multiple lumen hub with minimal volume of common fluid space

- Drugs do not mix in the lines
- Reduces the risk of drug interaction

Integral Anti-Reflux valve

- Prevents back-tracking of fluids and drugs

Polyurethane gravity line with Anti-Reflux valve

- Robust, kink resistant construction prevents occlusions
- Provided to allow rapid infusions