

MEDFUSER® FastFlow Sizes and Flow Rates

MedFuser® FastFlow elastomeric pumps are available in the following configurations

Pump size (ml)	Infusion rate (ml/h)	Duration	Pump size (ml)	Infusion rate (ml/h)	Duration
60	50	1h 12min	275	50	5h 30min
	100	36min		100	2h 45min
	200	18min		200	1h 22min
	250	14.4min		250	1h 6min
100	50	2h	450	50	9h
	100	1h		100	4h 30min
	200	30min		200	2h 15min
	250	24min		250	1h 48min
180	50	3h 36min	600	50	12h
	100	1h 48min		100	6h
	200	54min		200	3h
	250	43min		250	2h 24min

Technical specifications for different sizes of MEDFUSER® pumps

Pump size (ml)	Minimum capacity (ml)	Nominal capacity (ml)	Maximum capacity (ml)	Total residual volume (ml)
60	36	60	70	1
100	60	100	110	2
180	120	180	210	2
275	165	275	310	3
450	270	450	480	3
600	360	600	650	3

NOTE: Specifics may differ depending on the pump model and configuration of accessories.