

MEDFUSER® ONCO Sizes and Flow Rates

MedFuser®
ONCO
pump sizes (ml):

- 60
- 100
- 180
- 275
- 450
- 600

Infusion rates available (ml/h):

- 0.5
- 1.0
- 1.5
- 2.0
- 2.5
- 3.0
- 4.0
- 5.0
- 6.0
- 7.0
- 8.0
- 10

Technical specifications for different sizes of MEDFUSER® pumps

Pump size (ml)	Minimum capacity (ml)	Nominal capacity (ml)	Maximum capacity (ml)	Total residual volume (ml)
60	36	60	70	1
100	60	100	110 (*130)	2
180	120	180	210	2
275	165	275	310	3
450	270	450	480	3
600	360	600	650	3

NOTE: Specifics may differ depending on the pump model and configuration of accessories.

***NOTE:** MedFuser® ONCO only