

MEDFUSER® PCA Sizes and Flow Rates

MedFuser®
PCA
pump sizes (ml):

MedFuser® PCA ambulatory elastomeric pumps for pain management are available in the following configurations:

Nominal Fill
Volume

Patient Bolus Size (ml)

Lock-out Time (minutes)



***NOTE:** MedFuser® pumps allowing the patient to receive an on-demand dose of pain medication, in addition to an underlying basal infusion are also available and are referred to as MedFuser® PCA+Bsl .

Technical specifications for different sizes of MEDFUSER® pumps

Pump size (ml)	Minimum capacity (ml)	Nominal capacity (ml)	Maximum capacity (ml)	Total residual volume (ml)
60	36	60	70	1
100	60	100	110	2
180	120	180	210	2
275	165	275	310	3
450	270	450	480	3
600	360	600	650	3

NOTE: Specifics may differ depending on the pump model and configuration of accessories.