

MEDFUSER® MultiFlow Sizes and Flow Rates

MedFuser®
MultiFlow
pump sizes (ml):

60

100

180

275

450

600

Available modules for adjustment of infusion rate:

8-Step module 0 - 3.5 ml/h

0

0.5

1

1.5

2

2.5

3

3.5

8-Step module 0 - 7 ml/h

0

1

2

3

4

5

6

7

8-Step module 0 - 14 ml/h

0

2

4

6

8

10

12

14

Technical specifications for different sizes of MEDFUSER® pumps

Pump size (ml)	Minimum capacity (ml)	Nominal capacity (ml)	Maximum capacity (ml)	Total residual volume (ml)
60	36	60	70	1
100	60	100	110 (*130)	2
180	120	180	210	2
275	165	275	310	3
450	270	450	480	3
600	360	600	650	3

NOTE: Specifics may differ depending on the pump model and configuration of accessories.

***NOTE:** MedFuser® ONCO only