

MEDFUSER® MultiFlow SPLIT Sizes and Flow Rates

MedFuser® MultiFlow SPLIT pumps are available in six sizes from 60 ml to 600 ml nominal fill volume. Each pump size is available with any of the below three MultiFlow control modules:

MedFuser®
MultiFlow SPLIT
pump sizes (ml):

Available modules for adjustment of infusion rate:

60
100
180
275
450
600



8-Step module 0 - 3.5 ml/h							
0	0.5	1	1.5	2	2.5	3	3.5
8-Step module 0 - 7 ml/h							
0	1	2	3	4	5	6	7
8-Step module 0 - 14 ml/h							
0	2	4	6	8	10	12	14

Each of two channels of MedFuser® MultiFlow SPLIT receives **HALF** of the selected infusion rate.

Technical specifications for different sizes of MEDFUSER® pumps

Pump size (ml)	Minimum capacity (ml)	Nominal capacity (ml)	Maximum capacity (ml)	Total residual volume (ml)
60	36	60	70	1
100	60	100	110	2
180	120	180	210	2
275	165	275	310	3
450	270	450	480	3
600	360	600	650	3

NOTE: Specifics may differ depending on the pump model and configuration of accessories.