

MEDFUSER® MultiFlow SPLIT Key Elements

Elastomeric (silicone) balloon
providing stable infusion pressure
and proven drug stability.


Accurate fluid level indicator
provides peace of mind for
clinicians and patients.

Prominent product information to help identify the appropriate pump for safe patient treatment.

Robust, UV light shielding protective case.

Luer-lock filling port inlet with cap.

Integrated
hanger and clip.



Rate adjustment key.

MultiFlow rate
adjustment module.

Removable slide clamps give clinician control on starting and stopping infusions.

Two distal luer-lock connections to patient line.

Distal-end protective cap automatically helps to prime the patient line without spills.

All MultiFlow pumps are provided with a set of two particle and air filters.

