


MEDFUSER® PCA MultiFlow Sizes and Flow Rates

MedFuser®
PCA MultiFlow
pump sizes (ml):

MedFuser® PCA MultiFlow ambulatory elastomeric pumps for pain management are available in the following configurations:



Technical specifications for different sizes of MEDFUSER® pumps

Pump size (ml)	Minimum capacity (ml)	Nominal capacity (ml)	Maximum capacity (ml)	Total residual volume (ml)
60	36	60	70	1
100	60	100	110	2
180	120	180	210	2
275	165	275	310	3
450	270	450	480	3
600	360	600	650	3

NOTE: Specifics may differ depending on the pump model and configuration of accessories.