

MEDFUSER® NRFit Sizes and Flow Rates

MedFuser® NRFit
pump sizes (ml):

60
100
180
275
450
600



All MedFuser® NRFit pumps are available
with the following fixed infusion rates (ml/h):

0.5	1	1.5	2	2.5	3	4	5	6	7
8	10	13	20	25	50	100	200	250	

Technical specifications for different sizes of MEDFUSER® pumps

Pump size (ml)	Minimum capacity (ml)	Nominal capacity (ml)	Maximum capacity (ml)	Total residual volume (ml)
60	36	60	70	1
100	60	100	110	2
180	120	180	210	2
275	165	275	310	3
450	270	450	480	3
600	360	600	650	3

NOTE: Specifics may differ depending on the pump model and configuration of accessories.