

MEDFUSER® DUO NRFit Sizes and Flow Rates

MedFuser® DUO NRFit pumps are available in six sizes from 60 ml to 600 ml nominal fill volume. Each pump size is available with any of the below fixed flow rate :

MedFuser®

DUO NRFit

pump sizes (ml):

Available fixed infusion rates (ml/h):

60
100
180
275
450
600



0.5	1	1.5	2	2.5	3	4
5	6	7	8	10	13	

Each of the two channels of MedFuser® DUO NRFit receives the **FULL** rate specified.

Technical specifications for different sizes of MEDFUSER® pumps

Pump size (ml)	Minimum capacity (ml)	Nominal capacity (ml)	Maximum capacity (ml)	Total residual volume (ml)
60	36	60	70	1
100	60	100	110	2
180	120	180	210	2
275	165	275	310	3
450	270	450	480	3
600	360	600	650	3

NOTE: Specifics may differ depending on the pump model and configuration of accessories.