

MEDFUSER® MultiFlow NRFit Sizes and Flow Rates

MedFuser®
MultiFlow NRFit
pump sizes (ml):

60
100
180
275
450
600



Available modules for adjustment of infusion rate:

8-Step module 0 - 3.5 ml/h							
0	0.5	1	1.5	2	2.5	3	3.5
8-Step module 0 - 7 ml/h							
0	1	2	3	4	5	6	7
8-Step module 0 - 14 ml/h							
0	2	4	6	8	10	12	14

Technical specifications for different sizes of MEDFUSER® pumps

Pump size (ml)	Minimum capacity (ml)		Nominal capacity (ml)		Maximum capacity (ml)		Total residual volume (ml)
60	36	...	60	...	70	...	1
100	60	...	100	...	110	...	2
180	120	...	180	...	210	...	2
275	165	...	275	...	310	...	3
450	270	...	450	...	480	...	3
600	360	...	600	...	650	...	3

NOTE: Specifics may differ depending on the pump model and configuration of accessories.