

MEDFUSER® MultiFlow SPLIT NRFit Sizes and Flow Rates

MedFuser® MultiFlow SPLIT NRFit pumps are available in six sizes from 60 ml to 600 ml nominal fill volume. Each pump size is available with any of the below three MultiFlow control modules:

MedFuser® MultiFlow
SPLIT NRFit
pump sizes (ml):

- 60
- 100
- 180
- 275
- 450
- 600



Available modules for adjustment of infusion rate:

8-Step module 0 - 3.5 ml/h

00.511.522.533.5

8-Step module 0 - 7 ml/h

01234567

8-Step module 0 - 14 ml/h

02468101214

Each of two channels of MedFuser® MultiFlow SPLIT receives **HALF** of the selected infusion rate.

Technical specifications for different sizes of MEDFUSER® pumps

Pump size (ml)	Minimum capacity (ml)	Nominal capacity (ml)	Maximum capacity (ml)	Total residual volume (ml)
60	36	60	70	1
100	60	100	110	2
180	120	180	210	2
275	165	275	310	3
450	270	450	480	3
600	360	600	650	3

NOTE: Specifics may differ depending on the pump model and configuration of accessories.