

MEDFUSER® PCA MultiFlow NRFit Key Elements

Elastomeric (silicone) balloon providing stable infusion pressure and proven drug stability.

Prominent product information to help identify the appropriate pump for safe patient treatment.

Accurate fluid level indicator provides peace of mind for clinicians and patients.

Luer-lock filling port inlet with cap.

Integrated hanger and clip.

All PCA pumps are provided with a set of two particle and air filters.

Removable slide clamp gives clinician control on starting and stopping infusions.

Rate adjustment key.

PCA module with MultiFlow rate control

Distal-end protective cap automatically helps to prime the patient line without spills.

Yellow colour connector and hydrophobic cap identify NRFit® connector, compatible only with regional anaesthesia devices.

